

Rotameters

Economical Field Rotameters



320 Series

- Easy-to-read scales
- Provisions for panel mounting
- Ball type
- Lightweight, rugged, and portable

Measuring Range (L/min)	Scale (Inches)	Accuracy*	Cat. No.	Qty.
0.05 to 0.5	2	5%	320-2A05	ea
0.1 to 1 [†]	4	3%	320-4A1	ea
0.4 to 5	4	3%	320-4A5	ea
2 to 20	4	3%	320-4A20L [#]	ea
4 to 50	4	3%	320-440	ea
3 to 30	4	3%	320-530	ea
3 to 30	4	3%	320-100 [‡]	ea

* Full scale

[†] Product scale is marked in CCM (cc/min), which is equivalent to ml/min.

[‡] Fittings adapted for QuickTake 30 sample pump

[#] Fittings adapted for use with BioLite

Pneumatic Test Kit

The Pneumatic Test Kit is an all-in-one kit for testing sample pump flow range compensation. The kit includes a rotameter that measures 1 to 5 L/min, a valve, and a magnetic gauge that reads from 0 to 80 inches water back pressure. Mounted in heavy-duty casing.



Cat. No. 224-6580

Standard Rotameters

- 100-mm scale
- Steel base plate/stand
- Float type
- NIST certificate available; contact SKC



393 Series

Measuring Range (L/min) ^Δ	Scale (mm)	Accuracy	Cat. No.	Qty.
0.02 to 0.25	100	2.5%	393-002025	ea
0.3 to 3.4	100	2.5%**	393-0334	ea
0.6 to 5	100	2.5%**	393-0650	ea
1 to 13	100	2.5%**	393-1130	ea
2 to 26	100	2.5%**	393-2260	ea

Accessories			
Calidaptor, for hands-free IOM calibration	391-01	ea	
Cowled Head Adapter	391-05	ea	

** Full scale-VDI/VDE Standard 3513/2:2008


^Δ Calibrated at 68 F (20 C) and 1 atmosphere

Need Annual Calibration?

SKC CAL will keep your calibration devices in optimum condition.

Learn more on pages 42-43!

Standard and Field Rotameter Accuracy vs. chek-mate Calibrator



For the chek-mate Calibrator see pages 38-39